OXYGEN/PRESSURE PORT CONNECTOR

If pressure readings and/or additional oxygen are required, an Oxygen/Pressure Port Connector is available (Code 900HC452).

1. Connect the large end of the Connector (T) to the Swivel (G) on the mask and the small end of the Connector to the main breathing tube (R).
2. Open the Port Cap(s) (P) and place the Pressure and/or Oxygen Tube (S) firmly over the Port(s) (Q).


(The two Ports are identical; either Port may be used for oxygen or pressure measurement. Only remove the Port Cap from Ports that are being used.)

- **Note:** At a fixed flow rate of supplemental oxygen flow, the inhaled oxygen concentration will vary, depending on the pressure settings, patient breathing pattern, mask selection, and the leak rate.

- **Warning:** If oxygen is used with CPAP/Bi-Level machine, the oxygen flow must be turned off when the CPAP/Bi-Level machine is not operating.

Explanation of Warning: When the CPAP/Bi-Level device is not in operation, and the oxygen flow is left on, oxygen delivered into the ventilator tubing may accumulate within the CPAP/Bi-Level machine enclosure. Oxygen accumulated in the CPAP/Bi-Level machine enclosure will create a risk of fire.


- **Warning:** Do not smoke while using or in the vicinity of someone using oxygen with this mask.

BIAS FLOW INFORMATION

Nasal masks for CPAP and Bi-Level have exhalation ports to expel the air that you exhale from the mask.

PRESSURE (cmH ₂ O)	3	5	7	9	11	13	15	17	19	21	23	25
FLOW (L/MIN)	19	25	30	34	38	42	45	48	51	54	56	59

RESISTANCE TO FLOW

Pressure drop through mask @ 50 l/min 0.1 cmH₂O
Pressure drop through mask @ 100 l/min 0.4 cmH₂O

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MULTI-PATIENT CLEANING INSTRUCTIONS

Your FlexiFit™ 407 Nasal Mask is provided patient ready. To reuse between patients please follow the guidelines below.

CLEANING/ DESOILING

- Before cleaning, remove visible soil by disassembling and thoroughly cleaning the mask with a small brush under running potable tap water.

FAILURE TO PROPERLY CLEAN THIS DEVICE MAY RESULT IN INADEQUATE DISINFECTION.

HIGH LEVEL DISINFECTION PARAMETERS

- NB: For autoclave cycles, components should be enclosed in a sterilization peel pouch.

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MASK FRAME	AUTOCLAVING Autoclave, 15 min exposure at 121°C (250°F).
	OR
	THERMAL DISINFECTION Soak parts in hot water bath for 30 minutes at 70°C (158°F).
	OR
SILICONE SEAL	STERRAD® 100S Follow manufacturer instructions for Sterrad® 100S Sterilization system.
	OR
ELBOW	CIDEX® OPA Follow manufacturer instructions for high level disinfection. Immerse for 12 minutes at 20°C (68°F). Rinse thoroughly in potable water.
	OR
SWIVEL	SPOROX II® Follow manufacturer instructions for high level disinfection. Immerse for 6 hours at 20°C (68°F). Rinse thoroughly in potable water.

HEADGEAR	<i>Hand wash the Headgear in pure soap dissolved in lukewarm water. Air dry out of direct sunlight</i>
----------	--

FOAM CUSHION	<i>Replace with new part.</i>
--------------	-------------------------------

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THE FLEXIFIT™ 407 NASAL MASK HAS BEEN VALIDATED TO WITHSTAND 20 CYCLES OF THE ABOVE PARAMETERS. DO NOT USE IF BADLY DAMAGED.

Disclaimer: It is the user's responsibility to qualify any deviations from the recommended methods of reprocessing.

CONTACT YOUR LOCAL HOMECARE DEALER FOR ADDITIONAL INFORMATION.

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Fisher & Paykel
HEALTHCARE

FLEXIFIT™ 407 PATIENT INSTRUCTIONS FOR USE

The FlexiFit™ 407 Nasal Mask is intended for multiple patient or single patient adult use by individuals who have been diagnosed by a physician as requiring CPAP or Bi-Level ventilator treatment.

FITTING YOUR MASK

1. Holding your mask in one hand, spread the Stretchgear™ Headgear (H) with your free hand just enough to fit comfortably over your head.
 2. Lightly hold the mask against your face covering the nose then slowly pull the Headgear over your head.
 3. Attach the Headgear to the Glider™ Strap (E).
 4. Pull the lower straps as low as possible down back of neck. Gently tighten straps to remove slack. Connect to CPAP/Bi-Level machine set at prescribed pressure.
- If necessary, tighten the lower straps to eliminate any leaks (Avoid overtightening the straps).
 - When you need to remove your mask, unhook the Stretchgear™ Headgear from the Glider™ Strap.

FITTING TIPS

- Forehead pads should rest gently on the forehead without top straps being pulled too tight.
- If there are any leaks after gently tightening the top and lower straps, pull the mask slightly out from the face while straps are still attached. This will allow the Silicone Seal to inflate with air. Put the mask back onto the face.

CLEANING YOUR MASK

To ensure continuing safety and comfort for your FlexiFit™ 407 Mask:

BEFORE FIRST USE:

- Hand wash the Stretchgear™ Headgear (H) in pure soap dissolved in lukewarm water. Do not soak for longer than 10 minutes then rinse in fresh water and allow to air dry out of direct sunlight.

BEFORE EACH USE:

- Inspect the mask for deterioration. Do not use if aged or damaged.

AFTER EACH USE:

1. Disassemble the Silicone Seal (C), Foam Cushion (D) and the Headgear (H) from the Mask Base (A).
2. Wash the Silicone Seal and Mask Base in pure soap dissolved in lukewarm water. Do not soak for longer than 10 minutes, then rinse in fresh water.
3. Wipe the Foam Cushion with a damp cloth. (Do not immerse the Foam Cushion in water.)
4. Leave all parts to dry out of direct sun light before reassembling.

EVERY 7 DAYS:

- Hand wash the Stretchgear™ Headgear (H) in pure soap dissolved in lukewarm water. Do not soak for longer than 10 minutes then rinse in fresh water and allow to air dry out of direct sunlight.

CAUTION

DO NOT WASH THE MASK IN THE DISHWASHER.
DO NOT CLEAN THE MASK WITH PRODUCTS CONTAINING ALCOHOL,
CHLORINE, ANTIBACTERIAL AGENTS OR BLEACH.
DO NOT STORE THE MASK IN DIRECT SUNLIGHT.
DO NOT SOAK THE MASK FOR LONGER THAN 10 MINUTES.

The above actions may deteriorate or damage the mask and shorten its life.
If your FlexiFit™ 407 weakens or cracks, seek replacement immediately.

ASSEMBLING YOUR MASK

Your FlexiFit™ 407 Mask is fully assembled and ready to use. However, when you need to clean your mask follow these instructions for reassembly.

1. Place Silicone Seal (C) face down onto flat surface. With thick side of Foam Cushion (D) facing down, place into Silicone Seal.
2. Starting with the slot at the apex of the Silicone Seal, assemble onto Mask Base (A). Push firmly together.
3. Lay the Stretchgear™ Headgear (H) flat with the loop side down. Position the assembled mask over the Headgear and attach the four straps into the corresponding slots in the Forehead Rest (F) and the Glider™ Strap (E). This can be done without undoing the Velcro™ tabs, by sliding the Headgear into the slots.

The FlexiFit™ 407 and accessories are Latex free.

Many people find CPAP therapy causes dryness of the nose and throat. Heated humidification can substantially reduce these symptoms.

WARNINGS

- This mask should be used only with CPAP or Bi-Level equipment recommended by your physician or respiratory therapist. A mask should not be used unless the CPAP or Bi-Level machine is turned on and operating properly.

Explanation of Warning: CPAP and Bi-Level machines are intended to be used with special masks which have exhalation ports to allow continuous flow of air out of the mask. When the CPAP/Bi-Level machine is turned on and functioning properly, new air from the CPAP/Bi-Level machine flushes the exhaled air out through the mask exhalation ports. However, when the CPAP/Bi-Level machine is not operating, enough fresh air will not be provided through the mask, and exhaled air may be rebreathed. Rebreathing of exhaled air for longer than several minutes can in some circumstances lead to suffocations.

- At low CPAP pressures the flow through the exhalation ports may be inadequate to clear all exhaled gas from the tubing. Some rebreathing may occur.
- Operating pressure range is 3-25 cm H₂O.
- Do not block the gas exhaust ports. Before use, always ensure that gas is flowing out the exhaust ports.
- Do not use the mask if vomiting or feeling nauseous.
- Discontinue using the mask if there is an allergic reaction to any part of the mask. Consult your physician if this occurs.
- If re-using the FlexiFit™ 407 between different patients, refer to the Multi-Patient Cleaning Instructions for the FlexiFit™ 407 Nasal Mask.

CAUTION

- Use the mask for its intended use only as directed in this booklet.
- US Federal Law restricts this device to sale by, or on order of a physician.
- Operating temperature range is 5-50°C (40-120°F).

US PATENT NUMBER 6 615 834
Other Patents Pending

FOR ADDITIONAL INFORMATION CONTACT YOUR LOCAL HOMECARE PROVIDER

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